

Code O.N.S



**Sterilization** By EtO gas

**Purpose** Built-in small Nebulizer with face mask enables Reliable, Efficient, and convenient Nebulization of Medication.

- Main Features**
- Soft, clear aerosol mask with anatomical form.
  - Wide, adjustable, elastic strap with sealed ends to preclude disassembly of strap.
  - Integrated nose clip for proper positioning of mask on the mouth and nasal area.
  - Small volume Nebulizer, with larger surface area provided for unique convex cone design, to create maximum capillary action.
  - Easy assembly / disassembly with screw-on cap provides a leak-proof seal.
  - Color-coded (Green) top fits either the aerosol mask, or a mouth T piece for oral Nebulizing.
  - Nebulizes 3cc within max.10 minutes in horizontal or vertical position.

**Packing**

Supplied sterile in individual peelable poly bag.