



Code(s) : 1/15

Main Features:

- * Tri-beveled & finely sharpened point.
- * Maximum comfort to patient & staff.
- * Consistent performance.
- * Standard thin size (28G x 23.5)
- * Light – angle point.
- * Excellent fit to most popular lancing devices.
- * Friendly colours.
- * Sterile.

Packing:

- * Supplied in small boxes of 25 - 50 - 100 & 200 pcs each.

Supplies on request:

- * Other needle gauge sizes are available on request.
- * Other types: Whole st. st, flat type, punched blood lancets are also available, peel packs Gamma sterilised, 200 pcs / box.