



Code(s) : Single 1/8 - Twin 1/9

Main Features :

- * Blood tapered spike with sharp piercing tip, and wide internal tract.
- * 20ml drip chamber including cylindrical blood filter conforming to DIN 58360
- * Precise roller clamp.
- * T shape injection site to assure safety to patient and staff.
- * 125 cm. transparent tubing.
- * Luer adaptor, or Threaded luer lock adaptor to secure the hermetic locking and stability of luer / needle hub.
- * Ultra thin-wall double beveled siliconized G19 x 1.5" needle on request.

Supplies on request :

- Blood transfusion sets with twin compartment drip chambers assuring easy priming, constant flow and avoiding air bubbles.
- Blood transfusion set integrated with IV set.

Packing :

- Economic HDPE pack assuring sterile internal tract.
- Standard blister pack (medical grade paper, PP / PE laminate) assuring sterility to both internal and external tract.

Sterilization : By EtO gas.