



Codes : Plain 1/67 - Cuffed 1/68

Main Features :

- * Thermosensitive tubing, easily adapted to respiratory tract.
- * Split-resistant radiopaque line.
- * Smooth, bull-nose tip ensures atraumatic intubation.
- * Smoothly connected inflation tube maximizes patient comfort. Soft, thin walled cuff with effective sealing.
- * I.D & O.D of tubing are marked on each set providing immediate identification.
- * High volume, Low pressure cuff.

Sizes : **Plain** **Cuffed**

Size Fr / Ch	10	12	14	16	22	24	26	28	30	32	34	36
I.D	2.5	3.0	3.5	4.0	5.5	6.0	6.5	7.0	7.5	8.0	8.5	9.0
O.D	3.3	4.0	4.7	5.3	7.3	8.0	8.7	9.3	10.0	10.7	11.3	12.0

Supplies on request :

- Other sizes not mentioned above.
- Oral curved, Nasal curved, and reinforced designs.

Packing :

- Individually packed in sterile pouches.
- Each 10 pcs in secondary duplex box.
- Each secondary boxes in Tertiary Corrugated Carton.

Sterilization: By EtO gas.

