

Extension Lines (Straight & Branched)

Code EL001, EL002, EL003

MOHReg 259/2004/1



Sterilization By EtO gas.

Main Features

- 3.0 X 4.1 mm. tubing, as per any requested length.
- Female luer lock adapter and cover to the container side.
- Male luer lock adapter to patient side.
- Injection site and/or on-off clamp on request.

Supplies on Request

- Extension line with Y connector providing two container lines and one patient line.
- Extension line with two Y connectors providing four container lines and one patient line.
- Any special design can be supplied at minimum supply time.

Packing

Standard blister pack (medical grade paper, PA/PE laminate) ensuring sterility to both internal and external tracts