



**Code(s) :** Venous 1/21 - Arterial 1/22 - Y 1/23

**Main Features :**

- \* Provided with Male Luer Lock at one end (patient side) and Female Luer Lock at the other end (container side) providing secure fitting to all standard equipment.
- \* 2 Protective caps to maintain the sterility of fluid path.
- \* Soft clear medical grade tubing:
  - Small bore (2.0 x 1.0 mm), High pressure for arterial line.
  - Large bore (4.1 x 3.0 mm), Low pressure for venous line.
- \* Tube length and volume as requested.

**Supplies on request :**

- Extension line with Y connector providing two container lines and one patient line.
- Extension line with two Y connectors providing four container lines and one patient line.
- Extension lines with clamps and /or injection sites.
- Extension lines with flow meter.

**Packing :**

- Economic HDPE pack assuring sterile internal tract.
- OR- Standard blister pack ( medical grade paper, PP / PE laminate ) assuring sterility to both internal and external tract.

**Sterilization :** By EtO gas.