



Code(s) :1/73

Function :

* For both internal feeding and medicine administration to premature babies and small infants.

Main Features :

- * Manufactured out of Medical-grade PVC with frosted surface to eliminate need for lubricant, and to be non-irritant to delicate mucosa.
- * Distal end with soft rounded tip for atraumatic introduction, 2 opposite, displaced eyes minimize risk of Dumping syndrome.
- * Proximal end with female Luer adaptor and cover for easy connection to feeding funnel or syringes, and hygienic protection against contamination.
- * Radio -opaque line throughout the tube full length.

Length Size: 38&52 cm.

Packing :

- Sterile, Individually straight packed in peel-packs, Then packs of 50 pcs.

Sterilization : By ETO gas.