


Catheterization Therapy

Feeding Tube

Code FC 52

MOHReg 92/2001/1



Sterilization By EtO gas

Sizes Length: 52 cm.

Purpose For both internal feeding and medicine administration to premature babies and small infants.

Varieties **Color Codes and Dimensions:**

Ch	4	5	6	8	10	12
O.D	1.3	1.7	2.0	2.7	3.3	4.0
Color	Red	Grey	Green	Blue	Black	White

Main Features

- Manufactured out of Medical -grade PVC with frosted surface to eliminate need for lubricant, and to be non-irritant to delicate mucosa.
- Distal **end** with **soft rounded tip** foratraumatic introduction, **2 opposite, displaced** eyes minimize risk **of** Dumping syndrome.
- Proximal **end** with **female Luer adaptor and cover for easy** c onnection to feeding funnel or syringes, and hygenic **protection against** contamination
- **Radio** -opaque line throughtout the **tube full length.**
- Marks & Numbers at **each cm** for clear **identification of position, start 5cm from**

the distal end.

Packing Sterile, Individually straight packed in peel-packs,
Secondary pcks of 50 pcs.