

**Code** FMIV001

**MOHReg** 437/2004/1



**Sterilization** By EtO gas.

**Main Features** Complete system with integrated I.V administration set and flow meter.  
( No need to use any extension lines)

Controller:

- Constant flow rate.
- Precise flow rate control.
- Easily adjustable, easily visible.
- Flow rate of 5 to 250 ml. per hour.
- Not affected with full bag or near-empty bag, bed height changes and /or patient position changes.

IV Set:

- Tapered 2 way spike with hydrophobic filter for air inlet.
- Soft clear drip chamber with 15 micron solution disk filter.
- Additional clamp for short term interruption of infusion.

- T shaped injection site to assure safety to patient and staff.
- Luer lock adapter to secure the hermetic locking and the stability of the luer/needle hub connection.
- Full length of 170cm assuring patient comfort.
- Ultra-thin wall, double bevelled ciliconised G21x1.5" needle.

### **Packing**

Standard blister pack (medical grade paper, PA/PE laminate) assuring sterility to both internal and external tracts.