



Code(s) : 1/1

Main Features :

- * Tapered spike with sharp piercing tip, Medical grade ABS.
- * Soft, clear drip chamber for easy priming.
- * Precise roller clamp.
- * Luer adaptor, or Threaded luer lock adaptor to secure the hermetic locking and stability of luer / needle hub.
- * Self-Sealing tubular injection site manufactured out of medical grade TPE. T or Y shaped injection site to assure safety to patient and staff.
- * 125 cm glass-clear tubing.
- * Ultra thin-wall double beveled siliconized G21 X 1.5" needle.

Supplies On-Request :

Infusion Sets with: solution disk filter 15 micron, 2-way spike for glass bottles, fine micro dropper 60 drops / ml, Y injection site, flash-ball, and with any tube length.

Special I.V Sets :

- Amber colored drip chamber and tube for the administration of light sensitive drugs.
- I.V sets suitable for infusion pumps.

Packing :

- Economic HDPE pack assuring sterile internal tract.
- OR - Standard blister pack (medical grade paper, PP/PE laminate) assuring sterility to both internal and external tracts.

Sterilization : By EtO gas.

