

**Code** IV200, IV300, (201, 202, 301, 302 & others)

**MOHReg** 5/2004/3



**Sterilization** By EtO gas.

**Main Features**

- Tapered spike with sharp piercing tip, Medical grade ABS.
- Soft, clear drip chamber for easy priming.
- Precise roller clamp.
- Luer adaptor, or Threaded luer lock adaptor to secure the hermetic locking and stability of luer / needle hub.
- Self-Sealing tubular injection site manufactured out of medical grade TPE, or T-shape injection site to assure safety to patient and staff.
- 125 cm glass-clear tubing.
- Ultra thin-wall double beveled siliconized G21 X 1.5" needle.

**Supplies on Request**

Infusion Sets with: solution disk filter 15 micron, 2-way spike for glass bottles, fine micro dropper 60 drops / ml, Y injection site, flash-ball, and with any tube length.

Special I.V Sets :

- Amber coloured drip chamber and tube for the administration of light sensitive drugs.
- I.V sets suitable for infusion

pumps.

**Packing**

- Economic HDPE pack assuring sterile internal tract.

OR

- Standard blister pack ( medical grade paper, PA/PE laminate )  
assuring sterility to both internal and external tracts.

---