

Respiratory Therapy

Nasal Canula

Code OM, G&P - NC, G&P

MOHReg 92/2001/1



Main Features

- Designed to provide maximum freedom to patient, by leaving his mouth free for nutritional and communicational purposes
- Twin-prong nasal tips are designed to ensure equal volume of Oxygen to : both air inlets.
- Two meter long clear PVC tube is attached to ensure continuous flow of oxygen.
- Proximal end is fitted with soft, clear, funnel-shaped connector that is easy to connect to any Oxygen source.

Packing Supplied Non-sterile, individual ly packed in poly bag.