



**Code(s) :** 1/66

**Main Features**

- Designed to provide maximum freedom to patient, by leaving his mouth free for nutritional and communicational purposes
- Twin-prong nasal tips are designed to ensure equal volume of Oxygen to : both air inlets.
- Two meter long clear PVC tube is attached to ensure continuous flow of oxygen.
- Proximal end is fitted with soft, clear, funnel-shaped connector that is easy to connect to any Oxygen source.

**Packing**

Supplied Non-sterile, individual ly packed in poly bag.

**Sterilization :** By EtO gas.