

**Nelaton Catheters  
and other simple PVC catheters**

**Code(s) :** 1/142

**Purpose of use :**

Intermittent bladder catheterization by introduction through the urethra.

**Main Features :**

- \* Frosted surface minimizes need for lubricant.
- \* Soft, rounded tip facilitates traumatic introduction.
- \* Lateral eyes compatible with the catheter lumen insure maximum drainage without being plugged.
- \* Color coded funnels for easy identification of Ch / Fr .

**Types :**

Nelaton, Tiemann and female catheter.

**Sizes :**

- Ch / Fr 6, 8, 10, 12, 14, 16, 18, 20, 22 and 24

**Other types of catheters :**

- \* Duodenal, levin and ryle tubes.
- \* Stomach and rectal tubes.
- \* Feeding tubes and umbilical catheters.
- \* Suction catheters and oxygen catheters .

**Packing :**

-Individually packed in special sterilization pouches.

**Sterilization :** By EtO gas.

