

## Catheterization Therapy

### Nelaton Catheters

**MOHReg** 92/2001/1



**Sterilization** By EtO gas.

**Sizes**

Ch / Fr 6, 8, 10, 12, 14, 16, 18, 20, 22 and 24.

**Related Products**

- Duodenal, levin and ryle tubes.
- Stomach and rectal tubes.
- Feeding tubes and umbilical catheters.
- Suction catheters and oxygen catheters .

**Purpose**

Intermittent bladder catheterization by introduction through the urethra.

**Main Features**

- Frosted surface minimizes need for lubricant.
- Soft, rounded tip facilitates traumatic introduction.
- Lateral eyes compatible with the catheter lumen insure maximum drainage without being plugged.
- Colour coded funnels for easy identification of Ch / Fr .

Types:  
Nelaton, Tiemann and female catheter.

**Packing**

Individually packed in special sterilization pouches