

**Code** O.V.S



**Sterilization** By EtO gas.

**Purpose** Single unit Venturi system that provides five accurate concentrations of oxygen (from 25 to 50 % ) " through 2 individual, color coded diluter jets.

**Main Features**

- Soft, transparent, aerosol mask with anatomical form. Ideal for long term use.
- Wide, adjustable, elastic strap with sealed ends to preclude disassembly of strap.
- Integrated nose clip for proper positioning of mask on the mouth and nasal area.
- Two color- coded diluter jets.
  - - Olive-green for 24,26,28 and 30% concentration.
  - - White for 35,40 and 50% concentration.
- Easy change of diluter jets in no time.
- Flexible, transparent, anti- c l ick tubing of two meter length.
- Both tube ends are fitted with soft funnel-shaped connector for easy connection to oxygen

source and diluterjets.

**Packing** Supplied sterile in individual  
peelable poly bag.