

**Oxygen Masks  
( For Convenient Oxygen Therapy )**

**Code(s) :** 1/58

**Main Features :**

**Mask:**

- \* Soft, Clear, Medical grade PVC mask, designed for convenient use.
- \* Provided with adjustable strap for proper positioning on the nasal area.
- \* Two meter long clear PVC tube is attached to ensure continuous flow of oxygen.
- \* Proximal end is fitted with soft, clear, funnel-shaped connector that is easy to connect to any Oxygen source.
- \* Available in Two sizes, for Adults and Children.

**Packing :**

- Supplied Non-sterile, individually packed in poly bag.

**Sterilization :** By EtO gas.