

## Respiratory Therapy

### Oxygen Mask

**Code** OM, G&P - NC, G&P



**Main Features**

- Soft, Clear, Medical grade PVC mask, designed for convenient use.
- Provided with adjustable strap for proper positioning on the nasal area.
- Two meter long clear PVC tube is attached to ensure continuous flow of oxygen,
- Proximal end is fitted with soft, clear, funnel-shaped connector that is easy to connect to any Oxygen source.
- Available in Two sizes, for Adults and Children.

**Packing** Supplied Non-sterile, individual ly packed in poly bag