

Dialysis Therapy

Peritoneal Dialysis Sets

Code PT001, PT002

MOHReg 262/2004/1



Sterilization By EtO gas.

Main Features The peritoneal set includes three ways:

- Way to container (single or twin) that consist of:
 1. Tapered spike(s) with sharp piercing tip(s).
 2. Soft clear drip chamber(s).
 3. Precise roller clamp(s).
- Way to patient that consist of:
 1. On-off clamp.
 2. Self sealing TPE tubing for waste dialysate sample drawing.
 3. Luer lock male adapter for catheter connection with cover.
- Way to drain that consist of:
 1. On-off clamp.
 2. Luer adapter with cover.

Supplies on Request Sets with long drip chambers, with solution filter and with any special

design can be supplied.

Packing

- Economic HDPE pack assuring sterile internal tract.
- Standard blister pack (medical grade paper, PA / PE laminate) assuring sterility to both internal and external tract.