



**Code(s) :** 1/26

**Main Features :**

- \* The pre-Operative scrub is a disposable sterile scrub brush / sponge for scrubbing of nails, hands and arms.
- \* It is supplied with a disposable nail cleaner that can also be used to puncture the soap pouch.
- \* It is available in plain, medi-soap, chlorohexidine gluconate, or povidone iodine soap which is contained in a pouch between the brush and the sponge.
- \* It prevents trauma to the skin around the nails, on hands, and arms.

**Packing :**

-Individually packed in special paper / laminate sterilization pouches.

**Sterilization :** By EtO gas.