

Catheterization Therapy

Rectal Tube

Code RT 40

MOHReg 92/2001/1



Sterilization By EtO gas

Sizes **Length Size:** 40 cm. **Sizes:**

<i>FG</i>	<i>22</i>	<i>25</i>	<i>30</i>	<i>32</i>	<i>35</i>
<i>Col. Code</i>	<i>Rose</i>	<i>Blue</i>	<i>Dark Grey</i>	<i>Brown</i>	<i>Light Grey</i>

Purpose For administration of medicine or draining of rectum by means of saline or oil for example.

- Main Features**
- Manufactured out of Medical-grade kink resistant PVC, with low friction surface.
 - Smooth, closed, round distal end and 2 large lateral eyes permit effective drainage.
 - frosted surface to reduce need for lubricant.
 - Color-coded funnels fits standard adaptors.

Packing Supplied sterile, Individually and straight packed in peelable soft sterilization pouches.