



Code(s) : 1/54

Function :

* For administration of medicine or draining of rectum by means of saline or oil for example.

Main Features :

- * Manufactured out of Medical-grade kink resistant PVC, with low friction surface.
- * Smooth, closed, round distal end and 2 large lateral eyes permit effective drainage.
- * frosted surface to reduce need for lubricant.
- * Color-coded funnels fits standard adaptors.

Length Size: 40 cm.

Sizes:

Packing :

- Supplied sterile, Individually and straight packed in peelable soft sterilization pouches.

Sterilization : By ETO gas.