



Codes: 1/40

Main Features :

- * The product consists of two clamps in between which runs a knife.
- * The clamp and cut function is performed by squeezing together with a three finger grip.
- * The tool will be separated to two parts:
 - A clamp (the white one) is left on the baby.
 - The cutter and clamp part (the green piece) is left with the placenta.
- * the procedure is carried out in a single operation and thus saving time.
- * Splash of blood from the umbilical cord is avoided, and thereby the danger of spreading of blood carried infection as HIV, Hepatitis B and C is reduced.
- * The clamp cut part remaining with the placenta can be opened at any time for blood sampling, study of the cut, and if necessary to drain the placenta.

Packing :

- Individually packed in special paper / laminate pouches.

Related product :

- Simple umbilical cord clamps are also available.

Sterilization: By EtO gas.

