

Code UCC & UC

MOHReg 26/2001/1



Sterilization By EtO gas.

Related Products Simple umbilical cord clamps are also available.

Main Features

- The product consists of two clamps in between which runs a knife.
- The clamp and cut function is performed by squeezing together with a three finger grip.
- The tool will be separated to two parts:
- A clamp (the white one) is left on the baby.
- The cutter and clamp part (the green piece) is left with the placenta.
- the procedure is carried out in a single operation and thus saving time.
- Splash of blood from the umbilical cord is avoided, and thereby the danger of spreading of blood carried infection as HIV, Hepatitis B and C is reduced.
- The clamp cut part remaining with the placenta can be opened at any time for blood sampling, study of the cut, and if necessary to drain the placenta.

Packing Individually packed in special paper/laminate pouches.