

Small Vein Set (Scalp Vein Set)

Code S.V0010.8.6.4

MOHReg 4/2005/1



Sterilization By EtO gas.

Main Features Our way to better I.V. Therapy. Light, easy to use, easy on patient, provides an ideal route for intermittent infusions. Also withdrawing multiple blood samples. Thanks to its miniature design and slim lines.

- Thin-will siliconised needles to provide larger pathway reducing the possibility of lumen occlusion. With sharp, lancet, bevel edge assuring easier penetration and reducing patient discomfort.
- Anti-slip wing with special locking points, to act as both built-in needle holder for easy vein-puncture, and as effective skin attachment during I.V. administration.
- Soft, flexible, non-toxic medical grade tubing capable of resisting kinks, curls and sets.
I.D varies from 2.5 to 1.0mm according to needle gauge.
- Soft PVC luer adaptor with built-in hermetic cover forms a leak-proof seal when desired, and provides an ideal route for intermittent infusion.

Colour code and specification:

Gauge	19G	21G	23G	25G	27G
Colour	Cream	Green	Sky Blue	Orange	Grey
O.D.	1.0	0.80	0.60	0.50	0.40

Supplies on

Request Other gauges can be arranged on request also any needle Length.

Packing

Standard blister pack (medical grade paper, PA/PE laminate) assuring sterility to both internal and external tracts.