

Code SW 001

MOHReg 264/2004/1



Sterilization By EtO gas.

Main Features

- Designed to precisely control the flow of solution from GLYCIN containers to the urethral catheter, resector, or cystoscope.
- Wide 4.8x6.8 tubing to assure high flow.
- Stout trochar assures good connection to the solution administration port.
- Connection to the catheter, resector, or cystoscope through soft, wide, flexible and translucent TPE tube.
- Step-wise multi-click clamp to regulate flow precisely.

Supplies on Request

Any special design can be supplied in short time.

Packing

- Economic HDPE pack assuring sterile internal tract.
OR
- Standard blister pack (medical grade paper, PA/PE laminate) assuring sterility to both internal and

external tracts.