

Three way Stopcock

Code SC 001, 002

MOHReg 436/2004/1



Sterilization By EtO gas.

Main Features

- **Fully** transparent **housing**.
- Continuous straight flow channels.
- Axially and radially movable **swivel lock for secure and easy connection**.
- Can be turned by 360°.
- **Pressure** tested, pneumatically **and hydrostatically**.
- Tight **covers** for all **ends** .
- Can be supplied alone or with 25cm. tubing (3x4.1 mm) ended with luer lock adaptor and cover .

Supplies on Request Any tube length for extension tubing, can be supplied.

Packing Standard blister pack (medical grade paper, PA/PE laminate) assuring sterility to both internal and external tracts.