

Blood Banking

Transfer Canula and Set

Code TC 001,002

MOHReg 438/2004/1



Sterilization By EtO gas.

Purpose Used for reconstruction of lyophilized drugs (specially plasma), injection of drug additives or mixing of two sterile fluids.

- Main Features**
- Two individual spikes with 20cm tube and clamp or single spike with two terminal ends.
 - Tapered spike with sharp piercing tip, medical grade ABS.
 - It allows secure fitting to all closure types.
 - It assures insignificant residual.

Packing Standard blister pack (medical grade paper, PA / PE laminate) assuring sterility to both internal and external tract.