



Code(s) : 1/57

Function :

* Designed for intermittent, or continuous access to the umbilical-artery of sick newborn or premature babies for monitoring blood pressure, medicine administration, or blood sampling.

Main Features :

- * Medical-grade PVC tubing with silk finish surface to facilitate the passage of catheter in the vessel.
- * Smooth round distal end helps easy cannulation
- * No lateral eyes, to eliminate the chance of clot formation in the blind space.
- * Proximal end is fitted with female Luer adaptor for convenient connection to I.V therapy equipment.
- * Marking at 5, 10, & 15cm. from the distal end helps to ascertain depth of insertion.
- * Radio-opaque line is provided for position control.

Length Size: 35 cm.

Sizes: FG: 3.5, 5, 6 & 8.

Packing :

- Packed in peelable pouches, sterile, ready for use.
- Secondly pack : Box of 50.

Sterilization : By ETO gas.