

Catheterization Therapy

Umbilical Catheters

Code U cath

MOHReg 92/2001/1



Sterilization By EtO gas.

Sizes Length Size: 35 cm. Sizes: FG: 3.5,5,6&8.

Uses Designed for intermittent, or continuous access to the umbilical-artery of sick newborn or premature babies for monitoring blood pressure, medicine administration, or blood sampling.

Price 0.00

Main Features

- Medical-grade PVC tubing with silk finish surface to facilitate the passage of catheter in the vessel.
- Smooth round distal end helps easy cannulation.
- No lateral eyes, to eliminate the chance of clot formation in the blind space
- Proximal end is fitted with female Luer adaptor for convenient connection to I.V therapy equipment.
- Markings at 5, 10, & 15cm. from the distal end helps to ascertain depth of insertion.
- Radio-opaque line is provided for position

control.

Packing

- Packed in peelable pouches, sterile, ready for use.
- Secondary pack : Box of 50.