



Code (s) : 1/10

Main Features :

- * Blood tapered spike with sharp piercing tip, and wide internal tract.
- * 100 or 150 ml calibrated burette for blood administration, including:
 - Air filter 0.2 micron with integrated cover.
 - Self sealing injection site.
 - Shut-off valve .
 - Wide dropper for blood.
- * 20ml.drip chamber including cylindrical blood filter conforming to DIN 58360
- * Two Precise roller clamp.
- * Extra medication T injection site to assure safety to patient and staff.
- * 150 cm. glass clear tubing.
- * Luer adaptor, or Threaded luer lock adaptor to secure the hermetic locking and stability of luer / needle hub with related cover.

Packing :

Available in two types of packaging :

- Simple poly bag.
- Window - carton over packed with PP laminate.

Sterilization : By EtO gas.