

Volume Control Blood Transfusion Set

Code BVC 001 , 002, 003

MOHReg 10/2005/1



Sterilization By EtO gas.

Main Features

- Blood tapered spike with sharp piercing tip, and wide internal tract.
- 100 or 150 ml calibrated burette for blood administration, including.
 1. Air filter 0.2 micron with integrated cover.
 2. Self sealing injection site.
 3. Shut-off valve .
 4. Wide dropper for blood.
- 20ml.drip chamber including cylindrical blood filter conforming to DIN 58360
- Two Precise roller clamp.
- Extra medication T injection site to assure safety to patient and staff.
- 150 cm. glass clear tubing.
- Luer adaptor, or Threaded luer lock adaptor to secure the hermetic locking and stability of luer / needle hub with related cover.

Packing Available in two types of

packaging :

- Simple poly bag.
- Window - carton over packed with PP laminate.