



Code(s) : 1/4

Main Features :

- * Tapered spike with sharp piercing tip.
- * 100 or 150 ml calibrated burette for I.V administration, including:
 - Air filter 0.2 micron with integrated cover.
 - Self sealing injection site.
 - Shut-off valve (fixed or floating).
 - Stainless steel dropper 60 drops = 1 ml.
 - 15 micron disk solution filter.
- * Two Precise roller clamp.
- * Extra medication T injection site.
- * Luer adaptor, or Threaded luer lock adaptor to secure the hermetic locking and stability of luer / needle hub.

Supplies On-Request :

- Economic vergins having no cover for air filter, no shut-off valve, no solution filter.
- Sets with dropper 20 drops = 1 ml.

Packing :

- Economic HDPE pack assuring sterile internal tract.
- Standard blister pack (medical grade paper, PP/PE laminate)

Sterilization : By EtO gas.