

Code VC 001, 002, 003

MOHReg 260/2004/1



Sterilization By EtO gas.

Main Features

- Tapered spike with sharp piercing tip.
- 100 or 150 ml calibrated burette for I.V administration, including:
 1. Air filter 0.2 micron with integrated cover.
 2. Self sealing injection site.
 3. Shut-off valve (fixed or floating).
 4. Stainless steel dropper 60 drops = 1 ml.
 5. 15 micron disk solution filter.
- Two Precise roller clamp.
- Extra medication T injection site.
- Luer adaptor, or Threaded luer lock adaptor to secure the hermetic locking and stability of luer / needle hub.

Supplies on Request

- Economic vergins having no cover for air filter, no shut-off valve, no solution filter.
- Sets with dropper 20 drops = 1 ml.
- Latest in technology shut-off ball.

Packing

- Economic HDPE pack assuring sterile internal tract.
- Standard blister pack (medical grade paper, PA/PE laminate)
- Window - carton over packed with PP laminate.